HOW TO USE Zensation Acupressure Mat™

FREE DOWNLOADABLE VIDEO
www.zenathletics.com/workoutvids

Acupressure Mat is based on the principles of acupressure and Japanese Shiatsu massage. Acupressure Mat stimulates specific reflex points throughout the body, releases blocked energy, ease tense muscles and creates deep mental and physical relaxation. Benefits include increased level of energy, reduced inflammation and pain, improved sleep, metabolism and circulation.

Acupressure Mat has 8,820 acupressure points that provides equal pressure on the body. You simply lay on it, ensuring that your "trouble areas" are stimulated by the stimulation points—and let the pressure of your body against the mat do the rest! The mat's stimulation point stimulate blood circulation, relax tense muscles and induce "feel good" hormones (endorphins and oxytocin) that release feelings of wellbeing and calm. It has helped so many people with back pain and neck and shoulder stiffness, sleeping issues, and general fatigue even depression. Mat does not claim to be substitute or cure for any medical condition.

Soothe your tense muscle and feel good anytime, anywhere

Lie flat on the mat for 35 minutes before sleeping to calm body.

Direct skin contact is recommended. Take deep breaths and relax onto the mat.

Initiates great sleep therefore aids in relaxing stiff neck, shoulder pain and back ache.

Relaxing tension in the neck and face

Instant facelift and fresh glow. Lie down 5 minutes each side of the face. If uncomfortable, use thin cloth.

To relax neck, roll down towel for neck support and place under Acupressure Mat. Rest your neck directly on the mat.

Improve blood circulation

Step barefoot onto the Acupressure Mat to increase blood circulation in relieving fatigue and sore feet.

CLEANING INSTRUCTIONS:

Remove inside foam mat and hand wash outside cotton cover, hang to dry.

Zensation® Athletics

Toll Free Consumer Support • Numéro de téléphone sans frais 1-877-415-3425
Zensation Athletics, PürAthletics and PürEarth products are distributed by
TriMax Sports, Vancouver B.C. Canada V5X 4V6
Les produits Zensation Athletics, PürAthletics et PürEarth sont distribués par
email: service@zenathletics.com | 9 am to 5 pm (PST) Mon - Fri

*Acupressure Mat is not a substitute to the medical treatment but a supplemental device only.

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